

Product data sheet

LoFric Cath - Kit

1. Company data

1.1. Manufacturer:

Name Astra Tech AB
 Address P.O. Box 14, SE-431 21 MÖLNDAL
 Country SWEDEN
 Website www.astratech.com

1.2. Distributor:

Name Astra Tech A/S
 Address Husby Allé 19, 2630 Taastrup
 Phone no. 0045 4371 3377
 Fax no. 0045 4371 7865
 E-mail astratech.dk@astratech.com
 Website www.astratech.dk

Quality management systems:

Manufacturer: ISO 9001

Distributor: _____

2. Product

2.1 Trade name LoFric Cath-Kit

2.2. List of product range:

See "Product range" (if Excel file)

"Click here" (if web based)

2.3. More information about the product

See "More information" (if Excel file)

"Click here" (if web based)

3. Product components

(For sterile products including the item package/sterile barrier)

Cardboard: <input type="text" value="N"/>	Paper: <input type="text" value="Y"/>	Glass: <input type="text" value="N"/>
Metal: <input type="text" value="N"/>	Latex: <input type="text" value="N"/>	Cotton: <input type="text" value="N"/>
Other: <u>TiO2 as pigment</u>		Organic: <input type="text" value="N"/>

3.2. Plastics:

Types of plastic: EBA, PVC, PVP, PE, PP PVC: Virgin suspension, Medical grade
 If PVC, state plasticiser: DEHP % of PVC: 33

3.3. Heavy metals, chlorine, etc.:

Heavy metals: Types: _____

Chlorine bleached material: Y, ECF in item package Recycled materials:

Are substances from the Environmental Agency's list of undesirable substances included? :

State substances: Cyclohexanone

4. Packaging components

(For sterile products **excluding** the item package/sterile barrier)

Cardboard: Paper: Glass:
Metal: Latex: Cotton:
Other: Organic:

4.2. Plastics:

Types of plastic: PVC:
If PVC, state plasticiser: _____ % of PVC: _____

4.3. Heavy metals, chlorine, etc.:

Heavy metals: Types: _____
Chlorine bleached material: Recycled materials:

State any substances from the Environmental Agency's list of undesirable substances: _____
Note: By agreement with Sinerfa this point takes effect only from 1 June 2003!!

5. CE marking

5.1. CE mark: No: Yes: Class:

6. Environmental information

6.1. Does a recognised env. label exist for the product area? Yes: No:

6.2. Does the product carry an environmental label? No: Yes: Label: _____

6.3. Environmental management systems: (State recognised environmental management systems, if any)

Manufacturer: Distributor: _____

7. Shelf life

Shelf life after date of manufacture / date of sterilisation:

8. Recommended storage and handling

State recommendations

9. Recommendations for waste handling

Recommendation = X	Incineration	Landfill	Recycling	Recycling label
Product	x			
Item package	x			
Inner package			x	
Transport packaging			x	

of product /product cate

List of product range

Manufacturer product no.	Distributor product no.	EAN no. sales unit	Description	Size	Packaging unit	No. of packs per trsp carton	No. of cartons per pallet	Ste- rile
910800	910800	7392532015014	LOFRIC C-KIT NEL 08	8	1			*
910800	910800	7392532015021	LOFRIC C-KIT NEL 08	8	25	2	70	*
911000	911000	7392532015137	LOFRIC C-KIT NEL 10	10	1			*
911000	911000	7392532015144	LOFRIC C-KIT NEL 10	10	25	2	70	*
911200	911200	7392532015250	LOFRIC C-KIT NEL 12	12	1			*
911200	911200	7392532015267	LOFRIC C-KIT NEL 12	12	25	2	70	*
911400	911400	7392532015373	LOFRIC C-KIT NEL 14	14	1			*
911400	911400	7392532015380	LOFRIC C-KIT NEL 14	14	25	2	70	*
911600	911600	7392532015496	LOFRIC C-KIT NEL 16	16	1			*
911600	911600	7392532015502	LOFRIC C-KIT NEL 16	16	25	2	70	*
911800	911800	7392532015595	LOFRIC C-KIT NEL 18	18	1			*
911800	911800	7392532015601	LOFRIC C-KIT NEL 18	18	25	2	70	*
930600	930600	7392532016370	LOFRIC C-KIT PAED 06	6	1			*
930600	930600	7392532016387	LOFRIC C-KIT PAED 06	6	25	2	110	*
930800	930800	7392532016479	LOFRIC C-KIT PAED 08	8	1			*
930800	930800	7392532016486	LOFRIC C-KIT PAED 08	8	25	2	110	*
931000	931000	7392532016578	LOFRIC C-KIT PAED 10	10	1			*
931000	931000	7392532016585	LOFRIC C-KIT PAED 10	10	25	2	110	*
950800	950800	7392532018077	LOFRIC C-KIT FEM 08	8	1			*
950800	950800	7392532018084	LOFRIC C-KIT FEM 08	8	25	2	110	*
951000	951000	7392532018176	LOFRIC C-KIT FEM 10	10	1			*
951000	951000	7392532018183	LOFRIC C-KIT FEM 10	10	25	2	110	*
951200	951200	7392532018275	LOFRIC C-KIT FEM 12	12	1			*

951200	951200	7392532018282	LOFRIC C-KIT FEM 12	12	25	2	110	*
951400	951400	7392532018374	LOFRIC C-KIT FEM 14	14	1			*
951400	951400	7392532018381	LOFRIC C-KIT FEM 14	14	25	2	110	*
971000	971000	7392532020070	LOFRIC CATH KIT TIEM CH10	10	1			*
971000	971000	7392532020087	LOFRIC CATH KIT TIEM CH10	10	25	2	70	*
971200	971200	7392532020216	LOFRIC CATH KIT TIEM CH12	12	1			*
971200	971200	7392532020223	LOFRIC CATH KIT TIEM CH12	12	25	2	70	*
971400	971400	7392532020353	LOFRIC CATH KIT TIEM CH14	14	1			*
971400	971400	7392532020360	LOFRIC CATH KIT TIEM CH14	14	25	2	70	*
971600	971600	7392532020476	LOFRIC CATH-CIT TIEM CH16	16	1			*
971600	971600	7392532020483	LOFRIC CATH-CIT TIEM CH16	16	25	2	70	*
971800	971800	7392532020612	LOFRIC CATH-CIT TIEM CH18	18	1			*
971800	971800	7392532020629	LOFRIC CATH-CIT TIEM CH18	18	25	2	70	*