

# Product data sheet

## LoFric Dila-Cath

### 1. Company data

#### 1.1. Manufacturer:

Name Astra Tech AB  
 Address P.O. Box 14, SE-431 21 MÖLNDAL  
 Country SWEDEN  
 Website [www.astratech.com](http://www.astratech.com)

#### 1.2. Distributor:

Name Astra Tech A/S  
 Address Husby Allé 19, 2630 Taastrup  
 Phone no. 0045 4371 3377  
 Fax no. 0045 4371 7865  
 E-mail [astratech.dk@astratech.com](mailto:astratech.dk@astratech.com)  
 Website [www.astratech.dk](http://www.astratech.dk)

#### Quality management systems:

Manufacturer: ISO 9001

Distributor: \_\_\_\_\_

### 2. Product

2.1 Trade name LoFric Dila-Cath

#### 2.2. List of product range:

**See "Product range" (if Excel file)**

**"Click here" (if web based)**

#### 2.3. More information about the product

**See "More information" (if Excel file)**

**"Click here" (if web based)**

### 3. Product components

(For sterile products including the item package/sterile barrier)

Cardboard: <input type="checkbox"/> N	Paper: <input type="checkbox"/> Y	Glass: <input type="checkbox"/> N
Metal: <input type="checkbox"/> N	Latex: <input type="checkbox"/> N	Cotton: <input type="checkbox"/> N
Other: <u>TiO2 as pigment</u>		Organic: <input type="checkbox"/> N

#### 3.2. Plastics:

Types of plastic: PVC, PVP, PE, PP, Polyacrylic adhesive      PVC: Virgin suspension, Medical grade  
 If PVC, state plasticiser: DEHP      % of PVC: 33

#### 3.3. Heavy metals, chlorine, etc.:

Heavy metals:  N      Types: \_\_\_\_\_

Chlorine bleached material: Y, ECF in item package      Recycled materials:  N

Are substances from the Environmental Agency's list of undesirable substances included? :  Y

State substances: Cyclohexanone

#### 4. Packaging components

(For sterile products **excluding** the item package/sterile barrier)

Cardboard:

Paper:

Glass:

Metal:

Latex:

Cotton:

Organic:

Other: Bentonite clay

#### 4.2. Plastics:

Types of plastic: PET, PE

PVC:

If PVC, state plasticiser: \_\_\_\_\_

% of PVC: \_\_\_\_\_

#### 4.3. Heavy metals, chlorine, etc.:

Heavy metals:

Types: \_\_\_\_\_

Chlorine bleached material:

Recycled materials:

State any substances from the Environmental Agency's list of undesirable substances: \_\_\_\_\_

*Note: By agreement with Sinerfa this point takes effect only from 1 June 2003!!*

#### 5. CE marking

5.1. CE mark: No:

Yes:

Class: 1S

#### 6. Environmental information

6.1. Does a recognised env. label exist for the product area?

Yes:

No:

6.2. Does the product carry an environmental label?

No:

Yes:

Label: \_\_\_\_\_

6.3. Environmental management systems: (State recognised environmental management systems, if any)

Manufacturer: ISO 14001

Distributor: \_\_\_\_\_

#### 7. Shelf life

Shelf life after date of manufacture / date of sterilisation:

#### 8. Recommended storage and handling

State recommendations Dry and at room temperature

#### 9. Recommendations for waste handling

Recommendation = X	Incineration	Landfill	Recycling	Recycling label
Product	x			
Item package	x			
Inner package			x	
Transport packaging			x	





**MEDICO**  
INDUSTRIEN

**SINERFA**

Sygehusenes Indkøbs- og  
Erfaringsudvekslingsforum

**me of product /product categ**

**List of product range**

Manufacturer product no.	Distributor product no.	EAN no. sales unit	Description	Size	Packaging unit	No. of packs per trsp carton	No. of cartons per pallet	Ste- rile
891600	891600	7392532012792	DILA-CATH NELATON CH 16	16	1			*
891600	891600	7392532012808	DILA-CATH NELATON CH 16	16	25	2	91	*
891800	891800	7392532012815	DILA-CATH NELATON CH 18	18	1			*
891800	891800	7392532012822	DILA-CATH NELATON CH 18	18	25	2	91	*