

Product data sheet

LoFric INSTI - CATH

1. Company data

1.1. Manufacturer:

Name Astra Tech AB
 Address P.O. Box 14, SE-431 21 MÖLNDAL
 Country SWEDEN
 Website www.astratech.com

1.2. Distributor:

Name Astra Tech A/S
 Address Husby Allé 19, 2630 Taastrup
 Phone no. 0045 4371 3377
 Fax no. 0045 4371 7865
 E-mail astratech.dk@astratech.com
 Website www.astratech.dk

Quality management systems:

Manufacturer: ISO 9001

Distributor: _____

2. Product

2.1 Trade name LoFric Insti-Cath

2.2. List of product range:

See "Product range" (if Excel file)

"Click here" (if web based)

2.3. More information about the product

See "More information" (if Excel file)

"Click here" (if web based)

3. Product components

(For sterile products including the item package/sterile barrier)

Cardboard: <input type="checkbox" value="N"/>	Paper: <input type="checkbox" value="Y"/>	Glass: <input type="checkbox" value="N"/>
Metal: <input type="checkbox" value="N"/>	Latex: <input type="checkbox" value="N"/>	Cotton: <input type="checkbox" value="N"/>
Other: <u>TiO2 as pigment</u>		Organic: <input type="checkbox" value="N"/>

3.2. Plastics:

Types of plastic: PVC, PVP, PE, PP, Terluc (ABS), Polyacrylic adhesive PVC: Virgin suspension, Medic
 If PVC, state plasticiser: DEHP % of PVC: 33

3.3. Heavy metals, chlorine, etc.:

Heavy metals: Types: _____

Chlorine bleached material: Y, ECF in item package Recycled materials:

Are substances from the Environmental Agency's list of undesirable substances included? :

State substances: Cyclohexanone

4. Packaging components

(For sterile products **excluding** the item package/sterile barrier)

Cardboard:

Paper:

Glass:

Metal:

Latex:

Cotton:

Organic:

Other: Bentonite clay

4.2. Plastics:

Types of plastic: PET, PE

PVC:

If PVC, state plasticiser: _____

% of PVC: _____

4.3. Heavy metals, chlorine, etc.:

Heavy metals:

Types: _____

Chlorine bleached material:

Recycled materials:

State any substances from the Environmental Agency's list of undesirable substances: _____

Note: By agreement with Sinerfa this point takes effect only from 1 June 2003!!

5. CE marking

5.1. CE mark: No:

Yes:

Class: 1S

6. Environmental information

6.1. Does a recognised env. label exist for the product area?

Yes:

No:

6.2. Does the product carry an environmental label?

No:

Yes:

Label: _____

6.3. Environmental management systems: (State recognised environmental management systems, if any)

Manufacturer: ISO 14001

Distributor: _____

7. Shelf life

Shelf life after date of manufacture / date of sterilisation:

8. Recommended storage and handling

State recommendations Dry and at room temperature

9. Recommendations for waste handling

Recommendation = X	Incineration	Landfill	Recycling	Recycling label
Product	x			
Item package	x			
Inner package			x	
Transport packaging			x	

sal grade

me of product /product categ

List of product range

Manufacturer product no.	Distributor product no.	EAN no. sales unit	Description	Size	Packaging unit	No. of packs per trsp carton	No. of cartons per pallet	Sterile
800800	800800	7392532007576	LOFRIC I/CATH NEL 08	8	1			*
800800	800800	7392532007583	LOFRIC I/CATH NEL 08	8	25	4	63	*
801000	801000	7392532007651	LOFRIC I/CATH NEL 10	10	1			*
801000	801000	7392532007668	LOFRIC I/CATH NEL 10	10	25	4	63	*
801200	801200	7392532007736	LOFRIC I/CATH NEL 12	12	1			*
801200	801200	7392532007743	LOFRIC I/CATH NEL 12	12	25	4	63	*
801400	801400	7392532007811	LOFRIC I/CATH NEL 14	14	1			*
801400	801400	7392532007828	LOFRIC I/CATH NEL 14	14	25	4	63	*
810800	810800	7392532007897	LOFRIC INSTI-CATH FEM CH 08	8	1			*
810800	810800	7392532007903	LOFRIC INSTI-CATH FEM CH 08	8	25	4	99	*
811000	811000	7392532007934	LOFRIC INSTI-CATH FEM CH 10	10	1			*
811000	811000	7392532007941	LOFRIC INSTI-CATH FEM CH 10	10	25	4	99	*
811200	811200	7392532007972	LOFRIC INSTI-CATH FEM CH 12	12	1			*
811200	811200	7392532007989	LOFRIC INSTI-CATH FEM CH 12	12	25	4	99	*
811400	811400	7392532008016	LOFRIC INSTI-CATH FEM CH 14	14	1			*
811400	811400	7392532008023	LOFRIC INSTI-CATH FEM CH 14	14	25	4	99	*
812800	812800	7392532008078	LOFRIC INSTI-CATH PED CH 08	8	1			*
812800	812800	7392532008085	LOFRIC INSTI-CATH PED CH 08	8	25	4	99	*
813000	813000	7392532008115	LOFRIC INSTI-CATH PED CH 10	10	1			*
813000	813000	7392532008122	LOFRIC INSTI-CATH PED CH 10	10	25	4	99	*
821000	821000	7392532008153	LOFRIC INSTI-CATH TIEMANN CH10	10	1			*
821000	821000	7392532008160	LOFRIC INSTI-CATH TIEMANN CH10	10	25	4	63	*
821200	821200	7392532008191	LOFRIC INSTI-CATH TIEMANN CH12	12	1			*

821200	821200	7392532008207	LOFRIC INSTI-CATH TIEMANN CH12	12	25	4	63	*
821400	821400	7392532008238	LOFRIC INSTI-CATH TIEMANN CH14	14	1			*
821400	821400	7392532008245	LOFRIC INSTI-CATH TIEMANN CH14	14	25	4	63	*