

Product data sheet

LoFric Plus

1. Company data

1.1. Manufacturer:

Name Astra Tech AB
 Address P.O. Box 14, SE-431 21 MÖLNDAL
 Country SWEDEN
 Website www.astratech.com

1.2. Distributor:

Name Astra Tech A/S
 Address Husby Allé 19, 2630 Taastrup
 Phone no. 0045 4371 3377
 Fax no. 0045 4371 7865
 E-mail astratech.dk@astratech.dk
 Website www.astratech.dk

Quality management systems:

Manufacturer: ISO 9001

Distributor: _____

2. Product

2.1 Trade name LoFric Plus

2.2. List of product range:

See "Product range" (if Excel file)

"Click here" (if web based)

2.3. More information about the product

See "More information" (if Excel file)

"Click here" (if web based)

3. Product components

(For sterile products including the item package/sterile barrier)

Cardboard: Paper: Glass:

Metal: Latex: Cotton:

Other: TiO2 as pigment Organic:

3.2. Plastics:

Types of plastic: PEBA, PVP, PE, PP, Cyanoacrylic glue, Polyacrylic adhesive PVC:

If PVC, state plasticiser: _____ % of PVC: _____

3.3. Heavy metals, chlorine, etc.:

Heavy metals: Types: _____

Chlorine bleached material: Y, ECF in item package Recycled materials:

Are substances from the Environmental Agency's list of undesirable substances included? :

State substances: _____

4. Packaging components

(For sterile products **excluding** the item package/sterile barrier)

Cardboard:

Paper:

Glass:

Metal:

Latex:

Cotton:

Organic:

Other: Bentonite clay

4.2. Plastics:

Types of plastic: PET, PE, PvAE

PVC:

If PVC, state plasticiser: _____

% of PVC: _____

4.3. Heavy metals, chlorine, etc.:

Heavy metals:

Types: _____

Chlorine bleached material:

Recycled materials:

State any substances from the Environmental Agency's list of undesirable substances: _____

Note: By agreement with Sinerfa this point takes effect only from 1 June 2003!!

5. CE marking

5.1. CE mark: No:

Yes:

Class: 1S

6. Environmental information

6.1. Does a recognised env. label exist for the product area?

Yes:

No:

6.2. Does the product carry an environmental label?

No:

Yes:

Label: _____

6.3. Environmental management systems: (State recognised environmental management systems, if any)

Manufacturer: ISO 14001

Distributor: _____

7. Shelf life

Shelf life after date of manufacture / date of sterilisation:

8. Recommended storage and handling

State recommendations Dry and at room temperature

9. Recommendations for waste handling

Recommendation = X	Incineration	Landfill	Recycling	Recycling label
Product	x			
Item package	x			
Inner package			x	
Transport packaging			x	

of product /product cat

List of product range

Manufacturer product no.	Distributor product no.	EAN no. sales unit	Description	Size	Packaging unit	No. of packs per trsp carton	No. of cartons per pallet	Ste- rile
903800	903800	7392532014475	LOFR PLUS NELATON CH 08	8	1			*
903800	903800	7392532014482	LOFR PLUS NELATON CH 08	8	25	4	63	*
904000	904000	7392532014536	LOFR PLUS NELATON CH 10	10	1			*
904000	904000	7392532014543	LOFR PLUS NELATON CH 10	10	25	4	63	*
904200	904200	7392532014598	LOFR PLUS NELATON CH 12	12	1			*
904200	904200	7392532014604	LOFR PLUS NELATON CH 12	12	25	4	63	*
904400	904400	7392532014659	LOFR PLUS NELATON CH 14	14	1			*
904400	904400	7392532014666	LOFR PLUS NELATON CH 14	14	25	4	63	*
904600	904600	7392532014710	LOFR PLUS NELATON CH 16	16	1			*
904600	904600	7392532014727	LOFR PLUS NELATON CH 16	16	25	4	63	*
904800	904800	7392532014772	LOFR PLUS NELATON CH 18	18	1			*
904800	904800	7392532014789	LOFR PLUS NELATON CH 18	18	25	4	63	*
905000	905000	7392532014833	LOFR PLUS NELATON CH 20	20	1			*
905000	905000	7392532014840	LOFR PLUS NELATON CH 20	20	25	4	63	*
905200	905200	7392532014895	LOFR PLUS NELATON CH 22	22	1			*
905200	905200	7392532014901	LOFR PLUS NELATON CH 22	22	25	4	63	*
905400	905400	7392532014956	LOFR PLUS NELATON CH 24	24	1			*
905400	905400	7392532014963	LOFR PLUS NELATON CH 24	24	25	4	63	*
923600	923600	7392532016196	LOFR PLUS PED CH 06	6	1			*
923600	923600	7392532016202	LOFR PLUS PED CH 06	6	25	4	99	*
923800	923800	7392532016257	LOFR PLUS PED CH 08	8	1			*
923800	923800	7392532016264	LOFR PLUS PED CH 08	8	25	4	99	*
924000	924000	7392532016318	LOFR PLUS PED CH 10	10	1			*

924000	924000	7392532016325	LOFR PLUS PED CH 10	10	25	4	99	*
943800	943800	7392532017711	LOFR PLUS FEMALE CH 08	8	1			*
943800	943800	7392532017728	LOFR PLUS FEMALE CH 08	8	25	4	99	*
944000	944000	7392532017773	LOFR PLUS FEMALE CH 10	10	1			*
944000	944000	7392532017780	LOFR PLUS FEMALE CH 10	10	25	4	99	*
944200	944200	7392532017834	LOFR PLUS FEMALE CH 12	12	1			*
944200	944200	7392532017841	LOFR PLUS FEMALE CH 12	12	25	4	99	*
944400	944400	7392532017896	LOFR PLUS FEMALE CH 14	14	1			*
944400	944400	7392532017902	LOFR PLUS FEMALE CH 14	14	25	4	99	*
944600	944600	7392532017957	LOFR PLUS FEMALE CH 16	16	1			*
944600	944600	7392532017964	LOFR PLUS FEMALE CH 16	16	25	4	99	*
944800	944800	7392532018015	LOFR PLUS FEMALE CH 18	18	1			*
944800	944800	7392532018022	LOFR PLUS FEMALE CH 18	18	25	4	99	*
983800	983800	7392532026904	LOFR PLUS FEMALE 15CM CH 08	8	1			*
983800	983800	7392532026898	LOFR PLUS FEMALE 15CM CH 08	8	25	4	99	*
984000	984000	7392532026928	LOFR PLUS FEMALE 15CM CH 10	10	1			*
984000	984000	7392532026911	LOFR PLUS FEMALE 15CM CH 10	10	25	4	99	*
984200	984200	7392532026942	LOFR PLUS FEMALE 15CM CH 12	12	1			*
984200	984200	7392532026935	LOFR PLUS FEMALE 15CM CH 12	12	25	4	99	*
984400	984400	7392532026966	LOFR PLUS FEMALE 15CM CH 14	14	1			*
984400	984400	7392532026959	LOFR PLUS FEMALE 15CM CH 14	14	25	4	99	*
984600	984600	7392532026980	LOFR PLUS FEMALE 15CM CH 16	16	1			*
984600	984600	7392532026973	LOFR PLUS FEMALE 15CM CH 16	16	25	4	99	*
984800	984800	7392532027000	LOFR PLUS FEMALE 15CM CH 18	18	1			*
984800	984800	7392532026997	LOFR PLUS FEMALE 15CM CH 18	18	25	4	99	*
993600	993600	7392532021497	LOFR PLUS PED 30 CM CH 06	6	1			*
993600	993600	7392532021503	LOFR PLUS PED 30 CM CH 06	6	25	4	63	*
993800	993800	7392532021558	LOFR PLUS PED 30 CM CH 08	8	1			*
993800	993800	7392532021565	LOFR PLUS PED 30 CM CH 08	8	25	4	63	*
994000	994000	7392532021619	LOFR PLUS PED 30 CM CH 10	10	1			*
994000	994000	7392532021626	LOFR PLUS PED 30 CM CH 10	10	25	4	63	*