

Product data sheet

LoFric

1. Company data

1.1. Manufacturer:

Name Astra Tech AB
 Address P.O. Box 14, SE-431 21 MÖLNDAL
 Country SWEDEN
 Website www.astratech.com

1.2. Distributor:

Name Astra Tech A/S
 Address Husby Allé 19, 2630 Taastrup
 Phone no. 0045 4371 3377
 Fax no. 0045 4371 7865
 E-mail astratech.dk@astratech.com
 Website www.astratech.dk

Quality management systems:

Manufacturer: ISO 9001

Distributor: _____

2. Product

2.1 Trade name LoFric

2.2. List of product range:

See "Product range" (if Excel file)

"Click here" (if web based)

2.3. More information about the product

See "More information" (if Excel file)

"Click here" (if web based)

3. Product components

(For sterile products including the item package/sterile barrier)

Cardboard: N Paper: Y Glass: N

Metal: N Latex: N Cotton: N

Other: TiO2 as pigment Organic: N

3.2. Plastics:

Types of plastic: PVC, PVP, PE, PP, Polyacrylic adhesive PVC: Virgin suspension

If PVC, state plasticiser: DEHP % of PVC: 33

3.3. Heavy metals, chlorine, etc.:

Heavy metals: N Types: _____

Chlorine bleached material: Y, ECF in item package Recycled materials: N

Are substances from the Environmental Agency's list of undesirable substances included? : Y

State substances: Cyclohexanone

4. Packaging components

(For sterile products **excluding** the item package/sterile barrier)

Cardboard: Paper: Glass:
Metal: Latex: Cotton:
Other: Organic:

4.2. Plastics:

Types of plastic: PVC:
If PVC, state plasticiser: _____ % of PVC: _____

4.3. Heavy metals, chlorine, etc.:

Heavy metals: Types: _____
Chlorine bleached material: Recycled materials:

State any substances from the Environmental Agency's list of undesirable substances: _____
Note: By agreement with Sinerfa this point takes effect only from 1 June 2003!!

5. CE marking

5.1. CE mark: No: Yes: Class:

6. Environmental information

6.1. Does a recognised env. label exist for the product area? Yes: No:

6.2. Does the product carry an environmental label? No: Yes: Label: _____

6.3. Environmental management systems: (State recognised environmental management systems, if any)

Manufacturer: Distributor: _____

7. Shelf life

Shelf life after date of manufacture / date of sterilisation:

8. Recommended storage and handling

State recommendations

9. Recommendations for waste handling

Recommendation = X	Incineration	Landfill	Recycling	Recycling label
Product	x			
Item package	x			
Inner package			x	
Transport packaging			x	

of product /product cat

List of product range

Manufacturer product no.	Distributor product no.	EAN no. sales unit	Description	Size	Packaging unit	No. of packs per trsp carton	No. of cartons per pallet	Ste- rile
900800	900800	7392532012877	LOFRIC NELATON CH 08	8	1			*
900800	900800	7392532012884	LOFRIC NELATON CH 08	8	25	4	63	*
901000	901000	7392532013058	LOFRIC NELATON CH10	10	1			*
901000	901000	7392532013065	LOFRIC NELATON CH10	10	25	4	63	*
901200	901200	7392532013232	LOFRIC NELATON CH12	12	1			*
901200	901200	7392532013249	LOFRIC NELATON CH12	12	25	4	63	*
901400	901400	7392532013430	LOFRIC NELATON CH 14	14	1			*
901400	901400	7392532013447	LOFRIC NELATON CH 14	14	25	4	63	*
901600	901600	7392532013638	LOFRIC NELATON 16	16	1			*
901600	901600	7392532013645	LOFRIC NELATON 16	16	25	4	63	*
901800	901800	7392532013812	LOFRIC NELATON CH 18	18	1			*
901800	901800	7392532013829	LOFRIC NELATON CH 18	18	25	4	63	*
902000	902000	7392532013997	LOFRIC NELATON CH20	20	1			*
902000	902000	7392532014000	LOFRIC NELATON CH20	20	25	4	63	*
902200	902200	7392532014154	LOFRIC NELATON CH22	22	1			*
902200	902200	7392532014161	LOFRIC NELATON CH22	22	25	4	63	*
902400	902400	7392532014314	LOFRIC NELATON CH24	24	1			*
902400	902400	7392532014321	LOFRIC NELATON CH24	24	25	4	63	*
920600	920600	7392532015694	LOFRIC PEDIATRIC CH06	6	1			*
920600	920600	7392532015700	LOFRIC PEDIATRIC CH06	6	25	4	99	*
920800	920800	7392532015878	LOFRIC PEDIATRIC CH08	8	1			*
920800	920800	7392532015885	LOFRIC PEDIATRIC CH08	8	25	4	99	*
921000	921000	7392532016035	LOFRIC PEDIATRIC CH10	10	1			*

921000	921000	7392532016042	LOFRIC PEDIATRIC CH10	10	25	4	99	*
940800	940800	7392532016677	LOFRIC FEMALE CH08	8	1			*
940800	940800	7392532016684	LOFRIC FEMALE CH08	8	25	4	99	*
941000	941000	7392532016851	LOFRIC FEMALE CH10	10	1			*
941000	941000	7392532016868	LOFRIC FEMALE CH10	10	25	4	99	*
941200	941200	7392532017032	LOFRIC FEMALE CH12	12	1			*
941200	941200	7392532017049	LOFRIC FEMALE CH12	12	25	4	99	*
941400	941400	7392532017216	LOFRIC FEMALE CH14	14	1			*
941400	941400	7392532017223	LOFRIC FEMALE CH14	14	25	4	99	*
941600	941600	7392532017391	LOFRIC FEMALE CH16	16	1			*
941600	941600	7392532017407	LOFRIC FEMALE CH16	16	25	4	99	*
941800	941800	7392532017551	LOFRIC FEMALE CH18	18	1			*
941800	941800	7392532017568	LOFRIC FEMALE CH18	18	25	4	99	*
961000	961000	7392532018657	LOFRIC TIEMANN CH10	10	1			*
961000	961000	7392532018664	LOFRIC TIEMANN CH10	10	25	4	63	*
961200	961200	7392532018817	LOFRIC TIEMANN CH12	12	1			*
961200	961200	7392532018824	LOFRIC TIEMANN CH12	12	25	4	63	*
961400	961400	7392532018992	LOFRIC TIEMANN CH14	14	1			*
961400	961400	7392532019005	LOFRIC TIEMANN CH14	14	25	4	63	*
961600	961600	7392532019173	LOFRIC TIEMANN CH16	16	1			*
961600	961600	7392532019180	LOFRIC TIEMANN CH16	16	25	4	63	*
961800	961800	7392532019357	LOFRIC TIEMANN CH18	18	1			*
961800	961800	7392532019364	LOFRIC TIEMANN CH18	18	25	4	63	*
962000	962000	7392532019531	LOFRIC TIEMANN CH20	20	1			*
962000	962000	7392532019548	LOFRIC TIEMANN CH20	20	25	4	63	*
980800	980800	7392532020759	LOFRIC FEMALE 15 CM CH 08	8	1			*
980800	980800	7392532020766	LOFRIC FEMALE 15 CM CH 08	8	25	4	99	*
981000	981000	7392532020834	LOFRIC FEMALE 15 CM CH 10	10	1			*
981000	981000	7392532020841	LOFRIC FEMALE 15 CM CH 10	10	25	4	99	*
981200	981200	7392532020872	LOFRIC FEMALE 15 CM CH 12	12	1			*
981200	981200	7392532020889	LOFRIC FEMALE 15 CM CH 12	12	25	4	99	*
981400	981400	7392532020919	LOFRIC FEMALE 15 CM CH 14	14	1			*
981400	981400	7392532020926	LOFRIC FEMALE 15 CM CH 14	14	25	4	99	*
990600	990600	7392532020957	LOFRIC PED. 30 CM CH 06	6	1			*
990600	990600	7392532020964	LOFRIC PED. 30 CM CH 06	6	25	4	63	*
990800	990800	7392532021138	LOFRIC PED. 30 CM CH 08	8	1			*

990800	990800	7392532021145	LOFRIC PED. 30 CM CH 08	8	25	4	63	*
991000	991000	7392532021312	LOFRIC PED. 30 CM CH 10	10	1			*
991000	991000	7392532021329	LOFRIC PED. 30 CM CH 10	10	25	4	63	*